

Remarks

By the above amendment to the specification, the correct capitalization and generic terminology for the tradename "SuperScript" have been provided at page 29 (see the enclosed "SuperScript™ II Reverse Transcriptase" brochure, Invitrogen, 2002). A clerical error in claim 17 has also been corrected. Additionally, kit claims 24-26 have been deleted to focus the issues to those relating to the method claims. Since the amendments reduce the issues and place the application in condition for allowance or better condition for appeal, Applicant respectfully requests their entry.

Applicant acknowledges with appreciation the withdrawal of the objection to the drawing figures set forth in the previous Office Action.

Applicant also appreciates the Examiner's clarification that the objection to the specification for the use of the trademark "SuperScript" was based on the lack of generic terminology accompanying it. Accordingly, the generic terminology "reverse transcriptase" has been expressly added and the precise capitalization of the trademark has been corrected. Applicant therefore requests withdrawal of the objection to the specification.

In the final Office Action, the Examiner maintained the rejection of claims 1, 2, 4, 5, 7-14, 16-24, and 26 under 35 U.S.C. § 102 as being clearly anticipated by Linsley et al. (US 6,271,002). Although the Examiner withdrew the Section 102 rejection of claims 3, 6, 15, and 25 based on the Linsley et al. reference, the Examiner newly rejected these claims under 35 U.S.C. § 103 as being unpatentable over Linsley et al. in view of Gu et al.

(US 6,436,677). Upon entry of the above amendment, the rejections of claims 24-26 will be rendered moot. With respect to the claims remaining after entry of the amendment, the rejections are in error for reasons discussed below.

As previously pointed out by Applicant, the Linsley et al. reference fails to teach or suggest a method as defined in any of claims 1-23 in which a second strand of DNA is synthesized by contacting a thermostable DNA polymerase, selected from Bst DNA polymerase large fragment, Bst DNA polymerase native enzyme, Taq DNA polymerase, Pfu DNA polymerase, Tgo DNA polymerase, Phi29 DNA polymerase plus trehalose, T7 DNA polymerase plus trehalose, Klenow fragment of E. coli polymerase plus trehalose, and native E. coli DNA polymerase I plus trehalose, with a thermostable RNase H under conditions conducive to thermostable DNA polymerase activity. Citing column 15, lines 50-59, and Figure 1 of Linsley et al., the Examiner argued that "[s]ince the amplification of the double stranded DNA is subject to amplification via Taq polymerase, *arguably*, the first strand DNA is contacted with Taq polymerase to generate a second strand DNA (or its complement)" and that "Linsley et al. do disclose that the second strand synthesis occurs between 37 and 55 °C" (final Office Action, page 4, emphasis added). The Linsley et al. teachings cited by the Examiner, however, do not pertain to the synthesis of double stranded DNA, which is described in the reference at column 15, line 66, through column 16, line 49. Rather, the use of Taq polymerase is taught by Linsley et al. in reference to PCR amplification of double stranded DNA according to methods well known in the art (see, e.g., col. 18, lines 62-64). Moreover, the use of a reaction temperature of between

60 and 90 °C, followed by cooling to a temperature between 37 and 55 °C, is described at col. 15, lines 51-56, in reference to the synthesis of first strand, not second strand, DNA.

In reference to the synthesis of second strand DNA, the Linsley et al. patent merely describes the use of conventional DNA polymerases such as E. coli DNA polymerase I, Klenow fragment of E. coli DNA polymerase I, and T4 DNA polymerase (see Linsley et al., e.g., col. 16, lines 11-14; col. 41, lines 27-29) and an incubation temperature such as 37°C (col. 16, lines 5-8).

Thus, from a proper reading of the Linsley et al. reference there can be no reasonable argument that it teaches or even suggests the presently claimed methods.

The teachings of Gu et al. fail to cure the deficiencies of the primary reference.

In rejecting claims 3, 6, 15, and 25 for obviousness, the Examiner cited the secondary reference for its disclosure of Bst DNA polymerase as a thermostable polymerase having reverse transcriptase activity. Since the primary reference fails to teach or suggest the use of a thermostable polymerase in synthesizing second strand DNA, there would have been no motivation for a person of ordinary skill in the art to have substituted the Gu et al. polymerase for a conventional DNA polymerase in Linsley et al.'s second strand synthesis and further to have elevated the incubation temperature for the second strand synthesis reaction. Consequently, the claims are patentable over the combined teachings of the Linsley et al. and Gu et al. references.


In view of the foregoing, Applicant requests allowance of the application.

Applicant also requests the Examiner to return an initialed Form PTO-1449 to confirm consideration of the reference submitted with the Supplemental Information Disclosure Statement dated November 24, 2003.

In the event any fees are required for the filing of this Amendment, please charge all necessary fees to Deposit Account No. 10-0750.

Respectfully submitted,

Date: January 6, 2004


Linda S. Evans
Reg. No. 33,873

Johnson & Johnson
One Johnson & Johnson Plaza
New Brunswick, New Jersey 08933-7003
(858) 320-3406



High performance RT for reliability in every experiment

To achieve the best results in your reverse transcriptase applications, you need the best reverse transcriptase. SuperScript™ II reverse transcriptase delivers the highest cDNA yields, the most full-length, and the greatest sensitivity of any RT enzyme. Combine SuperScript™ II RT with your choice of optimized buffers, protocols, and kits to get high performance and consistent reliability in a range of applications.

High yields of full-length cDNA

SuperScript™ II reverse transcriptase⁵ is a point mutant of Moloney murine leukemia virus RT (M-MLV RT) engineered to be RNase H⁻. Like M-MLV, it is a DNA polymerase that synthesizes a complementary DNA strand from single-

stranded DNA, RNA, or a RNA:DNA hybrid. But compared with wild type RTs, SuperScript™ II delivers higher yields of full-length cDNA—superior performance that enhances results in all your experiments.



Eliminating RNase H activity is the key

In wild type RTs such as M-MLV and AMV, RNase H activity competes with 5'→3' DNA polymerase activity (Table 1) for the hybrid formed between the RNA template and the DNA primer. This competition inhibits the elongation of template-primer complex, resulting in lower yields of cDNA (Figure 1). RNase H activity also acts to diminish the yield and reduce the size of the cDNA by hydrolyzing

the RNA template as chain growth occurs. The elimination of RNase H activity allows higher yields and significantly greater full-length cDNA (1). You'll typically achieve > 50% more full-length cDNA synthesis and substantially greater yields of first-strand synthesis than obtained with AMV RT or M-MLV RT.


Table 1 - Reverse transcriptase activities

Reverse Transcriptase Activities		
Enzyme	DNA Polymerase	RNase H*
AMV RT	46,000	494
M-MLV RT	350,000	30
SuperScript™ II RT	320,000	not detected

Relative activities of RTs. Numbers represent units/mg. A unit is a nmole substrate incorporated or nmole substrate solubilized in 10 min. Both polymerase and RNase H activities listed were determined with MgCl₂ (as used in first-strand cDNA synthesis).

* For RNase H activity, the substrate used was [³H](A)n•(dT)n. For AMV RT, optimal unit reaction conditions use MgCl₂. For M-MLV RT, optimal unit reaction conditions use MnCl₂ (not MgCl₂). With MnCl₂, M-MLV RT has an RNase H activity of 2,670 and SuperScript™ II RT has an activity of 0.0018.

Figure 1 - More full-length, first-strand cDNA with SuperScript™ II



cDNA was synthesized in duplicate according to manufacturers' instructions from a mixture of 0.5 µg each of MAP-4 mRNA (5.2 kb), MAP-4/APC hybrid mRNA (9.5 kb), and 7.5 kb mRNA. Reactions were incubated at recommended temperatures.

Lanes 1: Supplier Q1
Lanes 2: Supplier Q2
Lanes 3: Supplier I
Lanes 4: SuperScript™ II


Higher thermostability for enhanced synthesis

The point mutations leading to inactivity of the RNase H domain in SuperScript™ II RT also improve the temperature stability, allowing first-strand synthesis at temperatures as high as 50°C (Figure 2).

Reverse transcription at higher temperatures improves

cDNA synthesis from RNA that has high GC content or secondary structure, or both. With higher thermostability than wild type RTs, SuperScript™ II lets you meet these challenges—while giving you the high yield and full length you need.

Figure 2 - Effect of incubation temperature on RT performance




An autoradiograph is shown of ³²P-labeled cDNA synthesized from a mixture of 0.25 µg each of 1.3-kb, 2.4-kb, 4.1-kb, and 7.5-kb RNA with 200 units of each RT at various temperatures.

Superior performance, even in long RT-PCR


The superior performance of SuperScript™ II extends to long RT-PCR. Although other RTs are able to synthesize full-length cDNA from long transcripts, SuperScript™ II RT outperforms them. First, more RNA is converted to cDNA, and second, significantly more of the cDNA synthesized is full length (Figure 3). The increase in the amount of full-length transcript generated can have a marked effect on the yield of long (> 5 kb) RT-PCR products (2,3). As Figure 4 shows, only SuperScript™ II RT-derived cDNA yielded a signal for the 6.8-kb transcript.

Figure 4 – Comparison of RTs for long RT-PCR



A comparison of RTs for long RT-PCR. Synthesis of cDNA was catalyzed by SuperScript™ II RT (S), M-MLV RT (M), or AMV RT (A). All samples were treated with RNase H. An amount equivalent to 0.5% of the starting RNA was amplified with Elongase® Enzyme Mix. (A) 5 µg total HeLa RNA for TSC 5.3 kb; (B) 1 µg total HeLa RNA for TSC 5.3 kb; (C) 5 µg total HeLa RNA for Pol ε 6.8 kb.

Figure 3 - Comparison of RTs for long RT-PCR



Comparison of first-strand cDNA synthesis. Each RT was used in the amount (10 units AMV RT and 200 units other RTs) and under conditions determined optimal for each enzyme using 1 µg of a 7.5-kb RNA template. Total yield of first-strand cDNA (■) was calculated from TCA precipitations analysis. Full-length synthesis (□) was determined following electrophoresis on an alkaline agarose gel and cutting and counting the full-length band and the rest of the lane. The counts of the full-length band were divided by total counts from the lane and multiplied by 100%.


Multiple kit formats, greater options

With SuperScript™ II RT, you have the ideal tool to prepare cDNA copies of mRNA for both polymerase chain reaction (PCR) or cDNA library construction (Figure 5). SuperScript™ II RT is available in the following formats:

- *Stand-alone enzyme* – reliable first-strand synthesis with flexibility for any protocol
- *First-strand synthesis kit* – pre-qualified, optimized reagents for the best first-strand synthesis, combined with flexibility in choosing a PCR enzyme
- *One-step RT-PCR kits* – convenient, single-tube RT-PCR for the highest sensitivity, suitable for higher throughput applications
- *Double-strand synthesis kit* – high-quality cDNA with flexibility for any protocol or cloning vector
- *Library construction kits* – complete kits for high-quality libraries

Whatever format you choose, you'll get the high yield and full length you need for all your downstream experiments.

Figure 5 – Process of first-strand synthesis to PCR or double-strand synthesis




Complete kit for high-quality cDNA, every time

The SuperScript™ First-Strand Synthesis System for RT-PCR^{4,5} is a complete kit optimized to generate high-quality first-strand cDNA from total or poly(A)⁺ RNA. Use the kit in conjunction with PCR to detect the presence of rare messages, quantify specific mRNA from small samples, or clone specific cDNA without cDNA library construction. With this format,


you can choose the PCR enzyme best suited for your application. You'll get consistently high first-strand yields and full-length cDNA up to 12 kb (Figure 6) as well as exceptional sensitivity. Detect as little as 50 to 100 molecules of RNA template from samples as small as 1.0 pg total RNA (Figure 7).

Figure 6 - Comparison of yield and length



cDNA was synthesized from 10 ng total HeLa RNA for CBP 1.6 kb, 100 ng total HeLa RNA for Pol ε 3.5 kb and 6.8 kb and TSC 5.3 kb, 0.5 µg total HeLa RNA for APC 8.9 kb, and 2 µg total rat brain RNA for Dynein 12.3 kb. One-tenth of the cDNA reaction was used for 35 cycles of PCR with Platinum® Taq DNA Polymerase High Fidelity. Lane M: λ DNA/Hind III Fragments.

Figure 7 - Sensitivity



cDNA was synthesized from 1.0, 10, 100, 1,000, and 10,000 pg total HeLa RNA using the SuperScript™ First-Strand Synthesis System for RT-PCR. One tenth of the cDNA reaction was amplified for 40 cycles with Platinum® Taq DNA Polymerase High Fidelity. Lane C: Control containing 10 ng of total HeLa RNA and no RT. M: Low DNA Mass™ Ladder.

One-step kits for unrivaled sensitivity

High performance. Convenient set-up. With SuperScript™ One-Step RT-PCR Systems ^{1,2,5,14}, you get both. These easy-to-use kits combine SuperScript™ II RT with Platinum® automatic hot-start PCR technology to deliver optimal RT-PCR results. Add total or poly(A)⁺ RNA to a pre-mixture of RT and PCR enzymes in a single tube, and you'll get the high yield, sensitivity, and specificity required for successful gene expression analysis. Two SuperScript™ One-Step RT-PCR Systems are available, both with streamlined protocols:


- *SuperScript™ One-Step RT-PCR System with Platinum® Taq DNA Polymerase.* For fast and easy screening of gene expression, detection of rare RNA and viral RNA, and

quantification of specific messages. You'll detect as little as 5 to 10 molecules of RNA template from as little as 0.1 pg of total RNA (Figure 8).

- *SuperScript™ One-Step RT-PCR System for Long Templates.* Made with Platinum® Taq DNA Polymerase High Fidelity, allowing you to amplify RNA targets up to 9 kb (Figure 9). It's the ideal kit for cloning and expression of your RT-PCR products.


With SuperScript™ One-Step Systems, you get a single-tube format to reduce reaction variability and lower risk of contamination. The improved reproducibility and rapid protocol facilitate higher throughput applications.

Figure 8 - Increased sensitivity with SuperScript™ One-Step RT-PCR System with Platinum® Taq DNA Polymerase.



RT-PCR reactions contained 0.01, 0.1, 1.0, 10, 100, and 1,000 pg total HeLa RNA (lanes 1-6, respectively) and gene specific primers to detect GAPDH 532 bp (Panel A) or β -actin 353 bp (Panel B). Lane C: Control. Lane M1: 1 Kb Plus DNA Ladder™. Lane M2: Low DNA Mass™ Ladder.

Figure 9 - Comparison of RT-PCR products, ranging in fragment size, from various RT-PCR systems.




Five different targets were amplified from 100 ng total HeLa RNA: Pol ϵ 338 bp, PP2A 1,093 bp, CBP 1,590 bp, RPA 2,057 bp, and Pol ϵ 2,980 bp (lanes 1-5 respectively). All reactions were run using the manufacturer's recommended reaction and cycling conditions. Lane M: 1 Kb Plus DNA Ladder™.

High-quality library construction

A winning strategy depends on a solid foundation. SuperScript™ II delivers the superior performance you need to construct high-quality cDNA libraries. With a choice of three formats, you're sure to find the optimal system for your needs:

- *SuperScript™ Choice System*⁵. Generate ready-to-clone double-stranded cDNA from a purified mRNA population (Figure 10). This complete kit contains everything you need, including a detailed protocol, to ensure you construct a high-quality library the first time.
- *SuperScript™ Plasmid Systems with Gateway™ Technology*^{5,19,45}. These systems include specialized vectors that enable you to construct Gateway™-compatible cDNA libraries with large numbers of full-length clones. For more information about the power of Gateway™ Technology, please visit our web site at www.invitrogen.com.
- *SuperScript™ Double-Stranded cDNA Synthesis Kit*⁵. Synthesize high-quality double-stranded cDNA from total or poly(A)⁺ RNA. Supply your RNA primer, adaptors, purification columns, cloning vector, and competent cells, for maximum flexibility in library design.

Figure 10 - Outline of procedure of SuperScript™ Choice System



Advanced RACE kit for full-length genes

The GeneRacer™ RLM RACE kit uses SuperScript™ II RT and an RNA ligase-mediated method to capture the full-length 5' end from long and complex mRNA. For more information, please visit our web site at www.invitrogen.com.

Use only the best RT

Why take needless risks with your samples? Capture the high performance and consistent reliability that only SuperScript™ II reverse transcriptase can deliver—and set the stage for success in all your downstream experiments.

Product	Quantity	Cat. no.
SuperScript™ II RT	2,000 U	18064-022
	10,000 U	18064-014
	4 x 10,000 U	18064-071
SuperScript™ First-Strand Synthesis System	50 rxns	11904-018
SuperScript™ One-Step RT-PCR System	25 rxns	10928-034
<i>with Platinum® Taq DNA Polymerase</i>	100 rxns	10928-042
SuperScript™ One-Step RT-PCR System for Long Templates	25 rxns	11922-010
<i>with Platinum® Taq DNA Polymerase High Fidelity</i>	100 rxns	11922-028
SuperScript™ Double-Stranded cDNA Synthesis Kit	10 rxns	11917-010
SuperScript™ Choice System for cDNA Synthesis	3 rxns	18090-019
SuperScript™ Plasmid System for cDNA Synthesis and Plasmid Cloning with Gateway™ Technology	3 rxns	18248-013
<i>with UltraMAX™ DH5α-FT™ Competent Cells</i>	1 kit	18248-047
<i>with MAX Efficiency® DH10B™ Competent Cells</i>	1 kit	19613-017
<i>with MAX Efficiency® DH5α™ Competent Cells</i>	1 kit	19617-018
<i>with ElectroMAX™ DH10B™ Competent Cells</i>	1 kit	19625-011
GeneRacer™ Kit		
<i>with SuperScript™ II RT and TOPO TA Cloning® Kit for Sequencing</i>	1 kit	L1502-01

References:

1. Gerard, G. *et al.* (1989) *Focus*® 11:66.
2. Westfall, B. *et al.* (1995) *Focus*® 17:62.
3. Nathan, M. *et al.* (1995) *Focus*® 17:78.

1,2,4,5,14,19,45 Products mentioned above are subject to the Limited Use Label License indicated by the superscript numbers. Please refer to the Invitrogen web site or catalog for the Limited Use Label License corresponding to the numbers indicated.



Printed in the U.S.A. ©2002 Invitrogen Corporation. Reproduction forbidden without permission.

Corporate headquarters:

1600 Faraday Avenue • Carlsbad, CA 92008 USA • Tel: 760 603 7200 • Fax: 760 602 6500 • Toll Free Tel: 800 955 6288 • E-mail: tech_service@invitrogen.com • www.invitrogen.com

European headquarters:

Invitrogen Ltd, 3 Fountain Drive • Inchinnan Business Park • Paisley PA4 9RF, UK • Tel: +44 (0) 41 814 6100 • Fax: +44 (0) 141 814 6260 • E-mail: eurotech@invitrogen.com